

CL-300

Computerized Lensmeter




Metro 
Ophthalmic

PERFORMANCE
YOU CAN COUNT ON

 **TOPCON**

Touch Screen Panel UV Transmittance Measurement Compact Slim Body

- Modern Design with Ease-of-Use -



CL-300 Computerized Lensmeter

With a new UV Measurement Function, the CL-300 provides comprehensive data to both eyeglass wearers and optometrists. Its compact, ergonomic design and intuitive touch screen panel ensures that measurements can be easily and quickly obtained.


Features

- » LCD Touch Panel
- » UV Transmittance Measurement
- » Compact Slim Body
- » Green Measurement Light Beam
- » Automatic Mono- & Multi-focal Detection

UV Transmittance Measurement Function

The CL-300's new UV Measurement Function provides information on the ultraviolet transmittance in the range of 0% to 100%, providing reliable measurement results to both eyeglass and sunglass wearers.

* UV Transmittance Measurement is performed with the light beam of 365 nm wavelength.



Multi-focal Lens Measurement

Not only single lenses but also multifocal lenses can be easily measured. The CL-300 also features automatic detection for single and progressive lenses.





example: Progressive Lens Measurement

The image shows a near point measurement screen. With simple operation, the ADD value can be obtained by following the guide instructed on this progressive measuring mode.

Green Light Beam


The CL-300 incorporates the green light reading beam, which is also known as e-line. Due to the wavelength of 546 nm, the e-line can be well-sensed by the human's eye and it is determined as a standard by ISO 7944. Using green light beam in our instrument enhances lens measurement precision and reduces the chance of error in the measurements.



CL-300 Features


Improved Lens Support

Due to its new design, the lens support of the CL-300 is 6 mm more in height compared to the conventional model. This change enables smooth and trouble-free measurement of high-curve lenses, preventing from troublesome contact of the lens edge with external parts of the lens support.




Large 5.7 inch LCD Monitor with Wider Directional Angle

The wide LCD monitor covers a 160 degree angle of observation allowing clear, compatible readings from virtually any position regardless of the operator's height.



White Ink Cartridge Standard

The newly featured white marking ink provides clear marking even on lenses with anti-staining and anti-fog coating. Due to these favorable qualities, TOPCON made the white marking ink standard in the new CL-300.



Models and Marking Ink

	White Ink	Red Ink
Cartridge Type	standard	optional
Holder Type	—	standard

* For availability in your region and other details please refer to the local TOPCON authorized representative.


Contact Lens Measurement

Contact lenses can be easily measured using the CL-300 lens support and ring accessories that come standard with the instrument.



Fast, Easy to Load Printer

Loading the thermal paper into the printer has never been easier. Just drop the printer paper into the receptacle and close the cover. The quick and quiet printer is ready to go. This improvement saves time and enhances workflow.



Printout Sample

TOPCON CL-300			
⟨R⟩	S	C	A
	+0.00	+0.00	180
ADD 0.75			
PSM	0.00	OUT	0.00
⟨L⟩	S	C	A
	+0.75	+0.25	85
ADD 0.75			
PSM	0.00	OUT	3.50 DWN
TOPCON Optical Store			

Smart & Smooth Operation


Easy and Fast Operation

The large, newly designed color LCD monitor with a 5.7-inch touch panel simplifies the interpretation of readings even with progressive and multi-focal lenses. The dynamic color measuring graphics follow the measuring process and present each step of the operation with different color, providing better control and ease of use.



1

Smooth lens table support enables easy placement of eyeglasses.



2

Right after the lens is placed, the optical center (+) appears with the "BLUE" bar.



3

The "ALIGNMENT OK" is displayed when the target (+) comes within the circle. When BEEP function is ON, the buzzer will sound.



4

The measured data is memorized automatically with the "MARKING OK" when the target is aligned in the exact center (at 0Δ).

In order to do marking, place down the lens retainer to secure the lens and marking is stamped accurately with the equipped marking cartridge.


Specifications

Measurement Range	Sphere	0 to ± 25 D (0.01/0.12/0.25D steps)
	Cylinder	0 to ± 10 D (0.01/0.12/0.25D steps)
	Axis	0 to 180°(1°step)
	Addition	0 to + 10D (0.01/0.12/0.25D steps)
	Prism	0 to 10 Δ
Cylinder Mode		MIX / - / +
Prism Mode		X-Y (Orthogonal Coordinates), P-B (Polar Coordinates), No Display
UV Transmittance		0 ~ 100% (1% / 5% steps)
Lens Detection		AUTO for single & progressive lenses
Contact Lens Measurement		Hard / Soft contact lens
Light Source	Refractive Measurement	Green LED 540 \pm 5nm
	UV Measurement	UV LED 365 \pm 5nm
Wavelength Criterion		e-line/d-line with ABBE compensation
Lens Diameter		ϕ 5 to 100 mm
Screen		5.7 inch color touch screen
PD Measurement Range*		50-90mm
External Interface		RS-232C , LAN *
Power Supply		100-240V / 50/60Hz 1.2A (AC Adapter) : AUTO power save
Dimensions		197(W) \times 220(D) \times 404(H) mm
Weight		3.8Kg

*CL-300PDL model only

Models	CL-300	CL-300PDL
PD measurement	—	included
Printer	included	included
LAN	—	included

System Chart



Contact Topcon subsidiaries or dealers for system configurations.

*Subject to change in design and/or specifications without advanced notice.

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.



TOPCON TOPCON CORPORATION

732-727-0607
www.metroophthalmic.com

Metro
Ophthalmic

PRINTED WITH SOY INK
© 2012-2013 TOPCON CORPORATION All Rights Reserved
Printed in Japan 2013 06-20 NP E070-3 TMS • TSM • AS