

OC-2200

Ophthalmic Chair



Signature
Series

KEY FEATURES:

- » Manual recline with a incorporated ambidextrous handle
- » Illuminated membrane switch control
- » Solid cast aluminum frame, armrest, footrest and base
- » Multi position full range articulated headrest
- » Available in a variety of upholstery colors
- » Single button folding armrest

PERFORMANCE
YOU CAN COUNT ON

Metro
Ophthalmic




TOPCON

Technical Specifications

Overall Dimensions	Width: 28"; Depth: 40"; Height: 55"
Elevation Stroke (Low Type)	9"
Elevation Stroke (High Type)	12.5"
Rotation	355 Degrees
Inclination Angle	70 Degrees
Footrest Height	3.5" / 1.5"
Full Recline Length	75"
Low Seat Height	Height: 20"; Stroke 9"; Top Dimension of 29"
High Seat Height	Height: 22.5" Stroke 12.5"; Top Dimension of 35"
Voltage & Frequency	120VA(±10%); 5A; 60Hz 240VA(±10%); 3A; 50/60Hz
Shipping Weight	281lb
Lift Capacity	500lb
ETL Listed	Complies with UL and CSA

Upholstery Selections

Standard Colors



Designer Series Colors



Signature
Series

