



520

Stylish Assisted Recline Tilt Chair

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

 **RELIANCE**
INTERNATIONAL

Metro 
Ophthalmic

520

Assisted Recline Tilt Exam Chair

Ergonomic Features

- Detached footswitch for hands-free raising/lowering.
- Large padded headrest keeps patients comfortable while reclined.
- Upholstered armrests rotate up to allow patients to move in and out from either chair side.
- Dual-rotation locks allow you to rotate chair 360° from a single position.
- Large hinged footrest safely supports patient's full weight, and folds away to save space.

Effortless Recline to 40°

- A manual recline bar, easily accessible from both sides of the chair, controls the assisted tilt of the chair up to 40°.

Power Hydraulic Base Boosts Efficiency

- The 520 chair's corded footswitch allows you to raise and lower the patient while maintaining contact.

Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Base available in #20 Pearl.
- All materials are fire retardant.
- Choose from 14 standard vinyl colors and endless custom selections – call Reliance at 1-800-735-0357 for a color card.

CSA Mark of Quality

- The 520 has passed all rigorous CSA, C/US tests, affirming its safety and durability and achieving certificate of compliance.

Tech Specs

- Low base height: 21 1/2" lowered, 32 1/2" raised. Base: 22 1/8" by 24 1/4". Depth upright: 43 1/2". Reclined: 59 1/2". Seat: 21 1/2" wide by 15" deep. Armrest distance: 21 5/8". Shipping weight: 310 lbs. Standard electrical: 120 volts, 60 Hz, 3.0 amperes. Optional: 230 volts, 50 Hz, 3.0 amperes.

