



980

Premium Full-Power Exam/Procedure Chair

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

 **RELIANCE**
INTERNATIONAL

Metro 
Ophthalmic

980

Premium Full-Power Exam/ Procedure Chair



Ergonomic Features

- Memory programming cuts wasted movement.
- Dual side controls eliminate the need to reach over a patient to access switches.
- Rotation locks allow you to swivel chair 330° from both sides.
- Armrests lock in place and release for easy patient access and security.
- Auto return button returns chair to upright and lowered position.
- Large hinged footrest supports patient's full weight for safe ingress and egress.
- Detached footswitch allows for hands-free patient positioning.
- Headrest adjusts with a single hand and locks into place. Pillow and alternating headrests also available.

Powered Hydraulic Base and Reclining Top

- Keep contact with patient with finger switches and a corded footswitch that raises, lowers or reclines the chair and auto return.

Long-Lasting Membrane Switches

- Backlit membrane switches, which last 10-times longer than conventional switches, contribute to the 980's durability. Eliminate need for high voltage wiring.
- The programmable safety function prevents children from operating chair.
- The memory programmable function automatically positions chair to preferred exam position.

Materials, Finishes and Upholstery

- Long lasting hand-sewn upholstery with high-density foam lasts for years.
- Removable cushions for easy cleaning or replacement.
- Selection of 14 standard colors and dozens of custom selections.
- Base available in #20 Pearl.
- Cast iron bases and hinges improve stability and durability.
- Call Reliance at 1-800-735-0357 for a color card.

The Mark of Quality

- The Model 980 has passed rigorous Laboratories testing, which assures quality and safety for you and your patients, achieving CSA, C/US certificate of compliance.

Tech Specs

- Diameter: 24 1/4", Depth of chair (without headrest) table position: 60", Seat: 18 3/4" wide, 18 1/4" deep. Two base heights. HIGH: 23 5/8" lowered, 34 1/4" raised. LOW 21 1/8" lowered, 29 5/8" raised. Distance between armrests: 19" Shipping Weight: 325 lbs. Standard electric: 120 volts, 50/60 Hz, 6.3 amperes total. Optional electric: 230 volts, 50 Hz, 4.0 amperes.